

Analysis of Six Aromatic Amines Using a Core Enhanced Technology Accucore HPLC Column

Instrument: [Jasco HPLC System](#)

Sample Preparation

Primary standards of each amine were prepared in acetonitrile at a concentration of 1 mg/mL, with the exception of 2,4,5 trimethylaniline which was supplied at a concentration of 200 ng/mL.

Working standard contained the following:

80 ng/mL of 2,4,5 trimethylaniline

200 µg/mL of 2,4-diaminotoluene, o-toluidine and 2-methoxy-5-methylaniline

60 µg/mL of 4,4-oxydianiline and 4,4-methylene-bis (2-chloroaniline)

Thermo Scientific Column

Part Number

Accucore PFP 2.6 µm 100 x 2.1mm

17426-102130

Measured pressure: 300 bar

Thermo Scientific Accela

Column temperature 45 °C

Injection volume 1 µL

Flow rate 0.6 mL/min

UV detection 254 nm

Mobile Phase


Mobile phase A: 25 mM Ammonium acetate pH 5.0

Mobile phase B: Acetonitrile

Gradient: 20-100%B in 1.5 minutes

Peak	Analyte	Retention time (min)
1	2,4-Diaminotoluene	0.616
2	4,4-Oxydianiline	1.048
3	O-Toluidine	1.105
4	2-Methoxy-5-methylaniline	1.214
5	2,4,5-Trimethylaniline	1.357
6	4,4-Methylene-bis (2-chloroaniline)	1.509

Table 1: Peak identity and position of the aromatic amines obtained from Accucore PFP



Reference: <https://fscimage.thermoscientific.com/images/D02480~.pdf>