

Analysis of Ylang Ylang Oil by GC-MS

Column: Zebron™ ZB-1MS, GC Cap. Column 10 m x 0.10 mm x 0.10 µm

Phase: 100% Dimethylpolysiloxane


Oven Profile: 60°C to 120°C at 15°C/min to 160°C at 5°C/min to 220°C at 20°C/min

Carrier Gas: Constant Flow Helium, 0.5 mL/min

Injection: Split 120:1, 0.2 µL @ 160°C

Detection: [Mass Selective \(MSD\) \(300°C\)](#) ([Ctrl + Click to follow link](#))

Sample Preparation Note: Sample is 10% oil in dichloromethane.



1. p-Methyl anisole
2. Methyl benzoate
3. Linalool
4. Benzyl acetate
5. Geraniol
6. Gerenyl acetate
7. Copaene
8. b-Caryophyllene
9. Cinnamyl acetate
10. Humulene
11. Germacrene
12. alpha-Amorphene
13. Farnesene
14. d-Cadinene
15. Farnesol
16. Benzyl benzoate
17. Benzyl salicylate
18. Farnesyl acetate

Reference: E. E. Stashenko, W. Torres and J. R. Martinez Morales, A study of the compositional variation of the essential oil of ylang ylang (*Cananga odorata* Hook. fil. et Thomson, forma genuina) during flower development. *J. High Resol. Chromatogr.*, 18, 101-104 (1995).